


CRISPR Library

Instructions



Use Instruction - Human Epigenetic Knockout Library

■ Product Info

This Human Epigenetic Knockout Library is applicable for research on human epigenetic-related genes, containing 20,051 gRNAs, targeting 2,508 genes. The vector backbone of the library is LentiCRISPRv2-Opti backbone which is the all-in-one vector system, that is, the gRNA and Cas9 gene are on the same vector.

■ Library Details

Product Name	Human Epigenetic Knockout Library
Product Catalog	YKO-Libr-H007
Product Details	<p>20,051 gRNAs (gRNA sequences see attachment);</p> <p>Single-plasmid system, can be directly used for library screening without first constructing Cas9 stable cell lines;</p> <p>Puromycin resistance, puromycin can be used for antibiotic screening upon cell infection;</p> <p>Plasmids paired with 3rd lentivirus packaging system can be directly used for virus packaging.</p> <p>*It is recommended to use Ubigen's Lentiviral Packaging Kit (Cat# YK-LVP-05)</p> <p>Targeting 2,508 genes, 8 gRNAs per gene;</p> <p>50 non-target control gRNA.</p>
Backbone Map	

Verification Primers	<p>LentiCRISPRv2-Opti-F: ATTTCTTGGGTAGTTTGCAGTTT</p> <p>LentiCRISPRv2-Opti-R: GACTCGGTGCCACTTTTTTCA</p>
Product Specifications	<p>Ready-to-use, endotoxin-free, maxiprep plasmids, verified by Next Generation Sequencing, with coverage >99% and uniformity <10.</p>

■ Product Use Instruction

Part 1. Lentivirus Packaging

Mix library plasmid constructs with 3rd generation lentiviral packaging constructs to be co-transfected into 293T cells (Recommend: Ubigene's 293T cell line specialized for virus packaging, cat#YC-A006). 48 or 72 hours upon transfection, collect lentiviral supernatant and the virus can be used upon concentration. The virus should be stored at -80°C.

Part 2. Library Plasmid Amplification

1. Library plasmid electroporation

Add 50 ng library plasmid to 25 µL electrocompetent cells with transformation efficiency $\geq 10^9$ cfu/ug, electroporate cells as per electroporation parameters. Upon electroporation, add 975 µL recovery medium, mix well and transfer to a tube, then add 1 ml recovery medium to the tube and mix well again. Repeat above steps one more time and get 2 electroporation end products, share the tubes at 37°C, 250 rpm

for 1 hour.

2. Culture of amplified library and calculation of transformation efficiency

1) Mix the 2 tubes of electroporation end products and take 10 μL and dilute with 990 μL recovery medium. Plate 20 μL dilution onto a 10 cm Petri dish and incubate plates at 32°C for 14 hours. Count the colonies in the dish. If the number of colonies is 20000X greater than 4×10^5 , move on to the next step. If it is less than 4×10^5 , redo this step.

2) For the remaining electroporation end product, evenly inoculate into 2 flasks each containing 500 mL LB+Amp liquid culture medium. Then incubate at 37°C for 14 hours.

3. Collect transformation products

1) Collect the bacteria to a 50 mL centrifuge tube

2) Centrifuge tubes to pellet bacteria, decant LB and weigh pellet (bacteria).

4. Maxiprep

Maxiprep the plasmid DNA according to the instruction for maxiprep kit, it is recommended to use endofree maxiprep kit from QIAGEN, MACHEREY-NAGEL, etc (e.g. EndoFree Plasmid Mega Kit from QIAGEN)

Part 3. Library Screen

1. Determine infect MOI

Dilute the library virus into different gradients, such as MOI=0.3, 0.5, 1, 5, 10, 30, 100 to infect the target cells (the cell confluency is 30-50%). Each gradient needs to be set with 2 wells. After 48 hours of infection, add puromycin according to the

settings in the table below for screening, and stop antibiotic screening when all cells in the blank group (cells not infected with virus) die. The MOI with a survival rate of 30% after antibiotic screening is the virus infection condition for the library screening experiment, that is, infect MOI (MOI=0.3 in some literatures actually refers to the amount of virus corresponding to the virus infection of 30% cells).

Group#	MOI	Antibiotic screening	Cell amount upon antibiotic screening	Survival rate upon antibiotic screening
Experimental group 1	0.3	Yes	N1	N1/M1
Experimental group 2	0.5	Yes	N2	N2/M2
Experimental group 3	1	Yes	N3	N3/M3
Experimental group 4	5	Yes	N4	N4/M4
Experimental group 5	10	Yes	N5	N5/M5
Experimental group 6	30	Yes	N6	N6/M6
Experimental group 7	100	Yes	N7	N7/M7
Infection blank group 1	0.3	No	M1	—
Infection blank group 2	0.5	No	M2	—
Infection blank group 3	1	No	M3	—
Infection blank group 4	5	No	M4	—
Infection blank group 5	10	No	M5	—
Infection blank group 6	30	No	M6	—
Infection blank group 7	100	No	M7	—
Blank group	0	Yes	—	—

2. Transduction of target cells with library virus

- ① Determine the amount of cells and virus

$$\text{Cell amount} = \frac{\text{gRNA\#} \times \text{gRNA coverage}}{30\%} \quad * \text{gRNA coverage} > 500 \text{ fold}$$

$$\text{Virus amount} = \text{cell amount} \times \text{infect MOI}$$

- ② Expand the cells according to the cell amount calculated in step ①, and prepare sufficient virus.

- ③ Use library virus infect the target cells, upon puromycin screening, divide the screened cells into experimental group and control group. Add target drugs to the experimental group for screening, upon screening, collect 1×10^7 cells and perform genome extraction for Next Generation Sequencing, and then compare and analyze the gRNAs of the experimental group and the control group.

■ Relevant products and service

Ubigen provides 35+ off-shelf libraries including Human/Mouse genome-wide plasmid library and some sub-libraries, and one-stop customized screening services for CRISPR-KO, CRISPRa, and CRISPRi including high-throughput sgRNA library construction, virus packaging, cell infection, drug screening, NGS sequencing, and data analysis, etc. Multiple deliverables fulfill different research needs!