(1) UBIGER

(1) UBIGENE **STR Authentication Report For Cell Line**

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us

Online Inquiry

1. Cell Name: AN3 CA

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer. JUBIGEN
- 3. Sample gene genotype analysis result:

| | Genotype | analysis resu | Its of STR and | Amelogenin l | oci in cells | |
|----------|--|---------------|----------------|--------------------------------|--------------|---------|
| STR Loci | STR profile for Sample cell Cell name: AN3 CA | | | STR profile for Cell Bank cell | | |
| | | | | Cell name: AN3-CA | | |
| | Allele1 | Allele2 | Allele3 | Allele1 | Allele2 | Allele3 |
| D5S818 | 11 | 14 | JUBIC | 11 | 14 | |
| D13S317 | 12 | 14 | | 12 | 14 | |
| D7S820 | 7.1 | 10 | | 7.1 | 10 | IBIGE |
| D16S539 | 10 | 14 | 15.0 | 10 | 14 | 15.0 |
| VWA | 14 | 20 | | 14 | 20 | |
| TH01 | 9.3 | 10 | | 9.3 | 10 | |
| AMEL | Х | Х | | Х | Х | |
| TPOX | 8 | 10 | .0 | 8 | 10 | |
| CSF1PO | 12 | 14 | 15.0 | 12 | 14 | 15.0 |
| D12S391 | 20 | 24 | | | | Œ, |
| FGA | 22 | 23 | | | | (1)BIO |
| D2S1338 | 20 | 23 | | | | |
| D21S11 | 29 | 30 | | | | |
| D18S51 | 15 | 18 | | | | |
| D8S1179 | 12 | 14 | | | | |

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us Online Inquiry

| D3S1358 | 17 | 17 | | | |
|---------|----|------|------|-----|----------|
| D6S1043 | 12 | 14 | 18.0 | | |
| PENTAE | 9 | 16 | | ENE | |
| D19S433 | 14 | 14 | UBIC | , | |
| PENTAD | 9 | 11 | | | -E1 |
| D1S1656 | 13 | 18.3 | | | * (1810) |

4. Conclusion: The STR results of <u>AN3 CA</u> cells are consistent with the genotypes of <u>AN3 CA</u> cell lines in <u>EXPASY</u> database, the cell ID corresponded to <u>CVCL_0028</u>, and the STR results <u>completely matched</u>. In the test, <u>Multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us
Online Inquiry

5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry