(1) UBIGER

## STR Authentication Report For Cell Line

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252



1. Cell Name: ARPE-19

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer. JUBIGEN
- 3. Sample gene genotype analysis result:

|          | Genotype  | e analysis resu | ılts of STR and | Amelogenin l                                      | oci in cells |         |  |
|----------|---|-----------------|-----------------|---|--------------|---------|--|
| STR Loci | STR profile for Sample cell  Cell name: ARPE-19 |                 |                 | STR profile for Cell Bank cell Cell name: ARPE-19 |              |         |  |
|          |   |                 |                 |   |              |         |  |
|          | D5S818  | 13              | 13              | NBIC  | 13           | 13      |  |
| D13S317  | 11  | 12              |                 | 11  | 12           |         |  |
| D7S820   | 9   | 11              |                 | 9   | 11           | - IBIGE |  |
| D16S539  | 9   | 11              |                 | 9   | 11           |         |  |
| VWA      | 16  | 19              |                 | 16  | 19           |         |  |
| TH01     | 6   | 9.3             |                 | 6   | 9.3          |         |  |
| AMEL     | Х   | Υ               |                 | Х   | Y            |         |  |
| TPOX     | 9   | 11              |                 | 9   | 11           |         |  |
| CSF1PO   | 11  | 11              | UBIC            | 11  | 11           |         |  |
| D12S391  | 21  | 22              |                 |   |              | a E     |  |
| FGA      | 23  | 23              |                 |   |              | NIBION  |  |
| D2S1338  | 19  | 19              |                 |   | (2)          |         |  |
| D21S11   | 28  | 29              |                 |   |              |         |  |
| D18S51   | 12  | 16              |                 |   |              |         |  |
| D8S1179  | 13  | 13              |                 |   |              |         |  |

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

| D3S1358 | 14 | 15 |      |     |       |
|---------|----|----|------|-----|-------|
| D6S1043 | 17 | 19 |      |     |       |
| PENTAE  | 7  | 11 | . C  | ENE |       |
| D19S433 | 12 | 13 | UBIC |     |       |
| PENTAD  | 11 | 13 |      |     | -E1   |
| D1S1656 | 11 | 15 |      |     | NIBIO |

4. Conclusion: The STR results of <u>ARPE-19</u> cells are consistent with the genotypes of <u>ARPE-19</u> cell lines in <u>DSMZ</u> database, the cell ID corresponded to <u>CRL-2302</u>, and the STR results <u>completely matched</u>. In the test, <u>No multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

## 5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry

3 UBIGENE