(1) UBIGER

STR Authentication Report For Cell Line

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252



1. Cell Name: BxPC-3

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
- 3. Sample gene genotype analysis result:

	Genotype	analysis resu	ılts of STR and	Amelogenin l	oci in cells		
STR Loci	STR profile for Sample cell Cell name: BxPC-3			STR profile for Cell Bank cell Cell name: BX-PC3			
	D5S818	11	11	BIC	11	11	
D13S317	11	11		11	11	-61	
D7S820	10	13		10	13	- IBIGE	
D16S539	9	11		9	11		
VWA	14	18		14	18		
TH01	9	9		9	9		
AMEL	Х	Х		Х	Х		
TPOX	8	8		8	8		
CSF1PO	13	13	UBIC	13	13		
D12S391	19.3	20				aE!	
FGA	20	21				NIBIO	
D2S1338	17	19					
D21S11	29	29					
D18S51	12	12					
D8S1179	13	13					

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

D3S1358	14	16			
D6S1043	12	12			
PENTAE	12	14		ENE	
D19S433	13	13	UBIC		
PENTAD	14	14			-E1
D1S1656	16	17.3			* I IBIOL

4. Conclusion: The STR results of <u>BxPC-3</u> cells are consistent with the genotypes of <u>BxPC-3</u> cell lines in <u>DSMZ</u> database, the cell ID corresponded to <u>760</u>, and the STR results <u>completely matched</u>. In the test, <u>No multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry

3 UBIGENE