(1) UBIGER

(1) UBIGENE **STR Authentication Report For Cell Line**

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us

Online Inquiry

1. Cell Name: C3A

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

	Genotype	analysis resu	Its of STR and	Amelogenin l	oci in cells		
STR Loci	STR profile for Sample cell Cell name: C3A			STR profile for Cell Bank cell Cell name: [HepG2/C3A, derivative of Hep G2 (ATCC HB-8065)]			
	D5S818	11	13		11	13	
D13S317	9	13		9	13	BIGH	
D7S820	10	10		10	10		
D16S539	12	13		12	13		
VWA	17	17		17	17		
TH01	9	9		9	9		
AMEL	Х	Y		X	Y		
TPOX	8	9	UBIC	8	9		
CSF1PO	10	11		10	11		
D12S391	21	25				- IBIGI	
FGA	22	25					
D2S1338	19	20					
D21S11	29	31					
D18S51	13	14					

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

D8S1179	15	16				
D3S1358	15	16				
D6S1043	13	13	. C	ENE		
PENTAE	15	20	UBIC			
D19S433	15.2	15.2				CEL
PENTAD	9	13				* IBIO
D1S1656	11	12			(2)	

4. Conclusion: The STR results of <u>C3A</u> cells are consistent with the genotypes of [HepG2/C3A, derivative of Hep G2 (ATCC HB-8065)] cell lines in <u>DSMZ</u> database, the cell ID corresponded to <u>CRL-10741</u>, and the STR results <u>completely matched</u>. In the test, **No multi allele was found** in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us
Online Inquiry