(1) UBIGER

## (1) UBIGENE **STR Authentication Report For Cell Line**

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us

Online Inquiry

1. Cell Name: CW-2

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer. JUBIGEN
- 3. Sample gene genotype analysis result:

	Genotype	analysis resu	ılts of STR and	Amelogenin l	oci in cells	
STR Loci	STR profile for Sample cell  Cell name: STR25012002			STR profile for Cell Bank cell Cell name: CW-2		
	D5S818	11	13	JUBIC	11	13
D13S317	9	11		9	11	
D7S820	10	13		10	13	IBIGE
D16S539	9	9		9	9	
VWA	18	21		18	21	
TH01	G 7	8		7	8	
AMEL	Х	Х		Х	Х	
TPOX	8	12		8	12	
CSF1PO	11	14	UBIC	11	14	
D12S391	18	19	20.0			Œ,
FGA	21	22				BIG
D2S1338	18	18				
D21S11	30	31				
D18S51	12	13				
D8S1179	15	17				

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us Online Inquiry

D3S1358	16	17			
D6S1043	16	16			
PENTAE	11	14	. C	ENE	
D19S433	12	13.2	14.2		
PENTAD	10	11			730
D1S1656	16	17	18.0		- IBIOL

4. Conclusion: The STR results of <u>CW-2</u> cells are consistent with the genotypes of <u>CW-2</u> cell lines in <u>EXPASY</u> database, the cell ID corresponded to <u>CVCL\_1151</u>, and the STR results <u>completely matched</u>. In the test, <u>No multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry

## 5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry