(1) UBIGER

## STR Authentication Report For Cell Line

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252



1. Cell Name: Caco-2

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
- 3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells								
STR Loci	STR profile for Sample cell  Cell name: Caco-2			STR profile for Cell Bank cell Cell name: CACO-2				
								Allele1
	D5S818	12	13	JUBIC	12	13		
D13S317	11	14		11	14	13.0		
D7S820	11	12		11	12	I IBIGE		
D16S539	12	12		12	13			
VWA	16	18		16	18			
TH01	6	6		6	6			
AMEL	Х	Х		Х	Х			
TPOX	9	11		9	11			
CSF1PO	11	11	OTUBIC	11	11			
D12S391	17	23				~E\		
FGA	19	19				BIGL		
D2S1338	17	25			(8)			
D21S11	30	30						
D18S51	12	12						
D8S1179	12	12						

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

D3S1358	G 17	17			
D6S1043	12	12			
PENTAE	7	7	. C	ENE	
D19S433	15	15	UBIC	,	
PENTAD	9	9			-E1
D1S1656	15	16			* IBIO

4. Conclusion: The STR results of <u>Caco-2</u> cells are consistent with the genotypes of <u>Caco-2</u> cell lines in <u>EXPASY</u> database, the cell ID corresponded to <u>CVCL\_0025</u>, and the STR results <u>completely matched</u>. In the test, <u>Multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

## 5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry

3 UBIGENE