(1) UBIGER

(1) UBIGENE **STR Authentication Report For Cell Line**

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us

Online Inquiry

1. Cell Name: EB-3

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
- 3. Sample gene genotype analysis result:

	Genotype	analysis resu	ults of STR and	Amelogenin l	oci in cells		
STR Loci	STR profile for Sample cell Cell name: STR24112711			STR profile for Cell Bank cell Cell name: EB-3 [EB3]			
	D5S818	9	10	JUBIC	9	10	
D13S317	12	14		12	14		
D7S820	11	11		11	11	1181GE	
D16S539	10	12		10	12		
VWA	17	19		17	19		
TH01	G 7	7		7	7		
AMEL	Х	Υ		Х	Υ		
ТРОХ	6	9	. C	6	9		
CSF1PO	12	15	NBIG	12	15		
D12S391	15	15					
FGA	22	22				I IBIGI	
D2S1338	17	22			(1)		
D21S11	29	29					
D18S51	15	17					
D8S1179	14	14					

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us Online Inquiry

D3S1358	15	16			
D6S1043	11	13			
PENTAE	14	16	. C	ENE	
D19S433	12.2	16.2	UBIC		
PENTAD	10	11			a E N
D1S1656	13	15			* IBIO

4. Conclusion: The STR results of <u>EB-3</u> cells are consistent with the genotypes of <u>EB-3</u> cell lines in <u>DSMZ</u> database, the cell ID corresponded to <u>CCL-85</u>, and the STR results <u>completely matched</u>. In the test, <u>No multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry



5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry

