(1) UBIGER

STR Authentication Report For Cell Line

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252





1. Cell Name: EL4

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 20- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
- 3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells									
	OGELAS	STR profile fo	or Sample cel	STR profile for Cell Bank cell					
STR Loci		Cell na	me: EL4	Cell name: EL4					
	Allele1	Allele2	Allele3	Allele4	Allele1	Allele2	Allele3		
18-3	18.0			1810	18.0				
4-2	20.3	21.3			20.3	21.3	-61		
6-7	15.0	16.0	17.0		15.0	16.0	17.0		
19-2	13.0	15.0			13.0	15.0			
1-2	15.0	18.0	20.0		15.0	18.0	20.0		
7-1	26.2	27.2			26.2	27.2			
8-1	16.0				16.0				
1-1	17.0	19.0		GEN	17.0	19.0			
3-2	14.0			JBIO	14.0				
2-1	16.0	17.0			16.0	17.0	a E.		
15-3	24.3	25.3	26.3		24.3	25.3	26.3		
6-4	18.0	19.0			18.0	19.0			
13-1	16.1	20.1	21.1		16.1	20.2	21.2		
11-2	16.0				16.0				
17-2	16.0				16.0				

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252



12-1	17.0				17.0		
5-5	17.0	18.0	20.0		17.0	18.0	20.0
X-1	28.0			GEN	28.0		
TH01				JBIO			
D5S818							~E\


Note: The cell lines were compared with the STR data of cell lines from ATCC, DSMZ, JCRB and RIKEN databases, the cell lines not included in the above cell banks could not be matched. D4S2408 and TH01 in the above sites is a human site, which is used to detect whether the cell is contaminated by human sources.

4. Conclusion: This cell line is identified as a mouse cell line. The STR results of <u>EL4</u> cells are consistent with the genotypes of <u>EL4</u> cell lines in <u>EXPASY</u> database, the cell ID corresponded to <u>CVCL_0255</u>, and the STR results <u>basically matched</u>. In the test, <u>Multiple</u> alleles were found, no cross-contamination, no human contamination in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252





Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

