(1) UBIGER

## STR Authentication Report For Cell Line

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252



1. Cell Name: HuH-7

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer. JUBIGEN
- 3. Sample gene genotype analysis result:

	Genotype	analysis resu	ılts of STR and	Amelogenin l	oci in cells		
STR Loci	STR profile for Sample cell  Cell name: HUH-7			STR profile for Cell Bank cell Cell name: HuH-7			
	D5S818	12	12	NBIC	12	12	
D13S317	10	11		10	11		
D7S820	11	11		11	11	- IBIGE	
D16S539	10	10		10	10		
VWA	16	18		16	18		
TH01	G 7	7		7	7		
AMEL	Х	Х		Х	Х		
TPOX	8	11		8	11		
CSF1PO	11	11	UBIC	11	11		
D12S391	20	21				CE!	
FGA	22	23				* I IBIO	
D2S1338	19	19					
D21S11	30	30					
D18S51	15	15					
D8S1179	14	14					

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

D3S1358	15	15			
D6S1043	13	15			
PENTAE	11	11	. C	ENE	
D19S433	13	14	UBIC		
PENTAD	12	12			-51
D1S1656	16	16			* I IBIO

4. Conclusion: The STR results of <u>HuH-7</u> cells are consistent with the genotypes of <u>HuH-7</u> cell lines in <u>DSMZ</u> database, the cell ID corresponded to <u>JCRB0403</u>, and the STR results <u>completely matched</u>. In the test, <u>No multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252



## 5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry

3 UBIGENE