(1) UBIGER

## (1) UBIGENE **STR Authentication Report For Cell Line**

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us

Online Inquiry

1. Cell Name: LN-229

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer. JUBIGEN
- 3. Sample gene genotype analysis result:

|          | Genotype   | analysis resu | ılts of STR and | Amelogenin l  | oci in cells |          |  |
|----------|--|---------------|-----------------|---|--------------|----------|--|
| STR Loci | STR profile for Sample cell  Cell name: LN-229-Y01-S |               |                 | STR profile for Cell Bank cell<br>Cell name: LN-229 |              |          |  |
|          |  |               |                 |   |              |          |  |
|          | D5S818   | 11            | 12              | NBIC  | 11           | 12       |  |
| D13S317  | 10   | 11            |                 | 10  | 11           |          |  |
| D7S820   | 8  | 11            |                 | 8   | 11           | - IBIGE  |  |
| D16S539  | 12   | 12            |                 | 12  | 12           |          |  |
| VWA      | 16   | 19            |                 | 16  | 19           |          |  |
| TH01     | 9.3  | 9.3           |                 | 9.3   | 9.3          |          |  |
| AMEL     | Х  | Х             |                 | Х   | Х            |          |  |
| TPOX     | 8  | 8             |                 | 8   | 8            |          |  |
| CSF1PO   | 12   | 12            | UBIC            | 12  | 12           |          |  |
| D12S391  | 20   | 22            |                 |   |              | Œ!       |  |
| FGA      | 23   | 23            |                 |   |              | * I IBIO |  |
| D2S1338  | 19   | 20            |                 |   | (2)          |          |  |
| D21S11   | 29   | 30            |                 |   |              |          |  |
| D18S51   | 13   | 15            |                 |   |              |          |  |
| D8S1179  | 13   | 14            |                 |   |              |          |  |

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us Online Inquiry

| D3S1358 | 16 | 17   |      |     |        |
|---------|----|------|------|-----|--------|
| D6S1043 | 11 | 11   |      |     |        |
| PENTAE  | 7  | 16   | . C  | ENE |        |
| D19S433 | 12 | 17.2 | UBIC |     |        |
| PENTAD  | 10 | 11   |      |     | -E1    |
| D1S1656 | 14 | 17   |      |     | * IBIO |

4. Conclusion: The STR results of <u>LN-229</u> cells are consistent with the genotypes of <u>LN-229</u> cell lines in <u>EXPASY</u> database, the cell ID corresponded to <u>CVCL\_0393</u>, and the STR results <u>completely matched</u>. In the test, <u>No multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us
Online Inquiry

## 5. Attached Image



Tel: +86 153 6067 3248(Int'I) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry