(1) UBIGER

(1) UBIGENE **STR Authentication Report For Cell Line**

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us

Online Inquiry

1. Cell Name: TTA1

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
- 3. Sample gene genotype analysis result:

	Genotype	analysis resu	ults of STR and	Amelogenin l	oci in cells		
STR Loci	STR profile for Sample cell Cell name: STR25030708			STR profile for Cell Bank cell Cell name: TTA-1			
	D5S818	12	13	NBIC	12	13	
D13S317	12	12		12	12	-61	
D7S820	11	12		11	12	- IBIGE	
D16S539	9	10		9	10		
VWA	17	17		17	17		
TH01	6	6		6	6		
AMEL	Х	Х		Х	Х		
TPOX	11	11		11	11		
CSF1PO	10	11	UBIC	10	11		
D12S391	18	18				a E	
FGA	23	23				* 1)BIO	
D2S1338	24	24			(2)		
D21S11	30	30					
D18S51	15	15					
D8S1179	13	15					

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us Online Inquiry

D3S1358	15	15			
D6S1043	14	14			
PENTAE	10	10	. C	ENE	
D19S433	14	14	UBIC		
PENTAD	13	13			~E\
D1S1656	16	16			* I IBIO

4. Conclusion: The STR results of <u>TTA1</u> cells are consistent with the genotypes of <u>TTA1</u> cell lines in <u>EXPASY</u> database, the cell ID corresponded to <u>CVCL_6297</u>, and the STR results <u>completely matched</u>. In the test, <u>No multi allele was found</u> in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us
Online Inquiry

5. Attached Image



Tel: +86 153 6067 3248(Int'I) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry