

STR Authentication Report For Cell Line

1. Cell Name: COLO-320

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: STR25120401			Cell name: COLO-320		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	12	12		12	12	
D13S317	11	11		11	11	
D7S820	9	12		9	12	
D16S539	11	12		11	12	
VWA	15	18		15	18	
TH01	9	9		9	9	
AMEL	X	X		X	X	
TPOX	8	9		8	9	
CSF1PO	11	11		11	11	
D12S391	18.3	19				
FGA	20	20				
D2S1338	18	25				
D21S11	33.2	33.2				
D18S51	15	15				
D8S1179	13	13				

D3S1358	17	17				
D6S1043	12	13				
PENTAE	11	11				
D19S433	16.2	16.2				
PENTAD	9	12				
D1S1656	16	17.3				

4. Conclusion: The STR results of **COLO-320** cells are consistent with the genotypes of **COLO-320** cell lines in **DSMZ** database, the cell ID corresponded to **ACC-144** , and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

5. Attached Image

