

STR Authentication Report For Cell Line

1. Cell Name: Li-7

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: Li-7-Y01-Z			Cell name: Li-7		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	9	11		9	11	
D13S317	11	11		11	11	
D7S820	10	11		10	11	
D16S539	9	9		9	9	
VWA	16	19		16	19	
TH01	7	7		7	7	
AMEL	X	X		X	X	
TPOX	8	11		8	11	
CSF1PO	12	12		12	12	
D12S391	22	24				
FGA	20	20				
D2S1338	17	19				
D21S11	30	30				
D18S51	16	16				
D8S1179	13	14				

D3S1358	15	16				
D6S1043	10	12				
PENTAE	5	5				
D19S433	12	14				
PENTAD	10	10				
D1S1656	13	17.3				

4. Conclusion: The STR results of **Li-7** cells are consistent with the genotypes of **Li-7** cell lines in **DSMZ** database, the cell ID corresponded to **RCB1941** , and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.