

STR Authentication Report For Cell Line

1. Cell Name: NCI-H23

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: NCI-H23-Y04-K			Cell name: NCI-H23 [H23]		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	12	13		12	13	
D13S317	12	12		12	12	
D7S820	9	10		9	10	
D16S539	11	11		11	11	
VWA	16	17		16	17	
TH01	6	6		6	6	
AMEL	X	X		X	X	
TPOX	8	9		8	9	
CSF1PO	10	10		10	10	
D12S391	15	17				
FGA	24	24				
D2S1338	18	23				
D21S11	30	30				
D18S51	14	14				
D8S1179	15	15				

D3S1358	15	15				
D6S1043	12	12				
PENTAE	7	17				
D19S433	12	14				
PENTAD	8	8				
D1S1656	14	14				

4. Conclusion: The STR results of **NCI-H23** cells are consistent with the genotypes of **NCI-H23** cell lines in **DSMZ** database, the cell ID corresponded to **CRL-5800** , and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

5. Attached Image

