

## STR Authentication Report For Cell Line

1. Cell Name: RAJI

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: Raji-Y01-Z			Cell name: RAJI'A'		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	10	13		10	13	
D13S317	13	13		13	13	
D7S820	10	10		10	10	
D16S539	8	11		8	11	
VWA	17	19		17	19	
TH01	6	7		6	7	
AMEL	X	Y		X	Y	
TPOX	8	13		8	13	
CSF1PO	10	12		10	12	
D12S391	18	19				
FGA	19	27				
D2S1338	22	22				
D21S11	28	31				
D18S51	17	17				
D8S1179	14	15				

D3S1358	15	16				
D6S1043	21	21				
PENTAE	5	13				
D19S433	14	14.2				
PENTAD	9	9				
D1S1656	14	14				

4. Conclusion: The STR results of **RAJI** cells are consistent with the genotypes of **RAJI** cell lines in **DSMZ** database, the cell ID corresponded to **HPACC** , and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

## 5. Attached Image

