

## STR Authentication Report For Cell Line

1. Cell Name: SU-DHL-6

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: SU-DHL-6			Cell name: SU-DHL-6		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	12	12		12	12	
D13S317	12	14		12	14	
D7S820	10	10		10	10	
D16S539	11	12		11	12	
VWA	14	17		14	17	
TH01	6	9.3		6	9.3	
AMEL	X	X		X	X	
TPOX	11	12		11	12	
CSF1PO	10	10		10	10	
D12S391	17	17				
FGA	20	24				
D2S1338	18	19				
D21S11	29	29				
D18S51	12	13				
D8S1179	14	14				

D3S1358	14	17				
D6S1043	13	13				
PENTAE	5	12				
D19S433	14	14				
PENTAD	8	10				
D1S1656	17.3	18.3				

4. Conclusion: The STR results of **SU-DHL-6** cells are consistent with the genotypes of **SU-DHL-6** cell lines in **DSMZ** database, the cell ID corresponded to **ACC-572** , and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

### 5. Attached Image

