

STR Authentication Report For Cell Line

1. Cell Name: SW620

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: SW620-Y01-FH			Cell name: SW620 [SW-620]		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	13	13		13	13	
D13S317	12	12		12	12	
D7S820	8	9		8	9	
D16S539	9	13		9	13	
VWA	16	16		16	16	
TH01	8	8		8	8	
AMEL	X	X		X	X	
TPOX	11	11		11	11	
CSF1PO	13	14		13	14	
D12S391	17	17				
FGA	24	24				
D2S1338	17	24				
D21S11	30	30.2				
D18S51	13	13				
D8S1179	13	13				

D3S1358	16	16				
D6S1043	11	12				
PENTAE	10	10				
D19S433	13	13				
PENTAD	6.2	9	15.0			
D1S1656	13	14				

4. Conclusion: The STR results of **SW620** cells are consistent with the genotypes of **SW620** cell lines in **DSMZ** database, the cell ID corresponded to **CCL-227** , and the STR results **completely matched**. In the test, **Multi allele was found** in this cell line.

5. Attached Image



