

## STR Authentication Report For Cell Line

1. Cell Name: TPC-1

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: TPC-1			Cell name: TPC-1		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	8	10		8	10	
D13S317	11	12		11	12	
D7S820	11	11		11	11	
D16S539	9	9		9	9	
VWA	14	18		14	18	
TH01	9	9.3		9	9	
AMEL	X	X		X	X	
TPOX	8	11		11	11	
CSF1PO	11	12		11	12	
D12S391	20	26				
FGA	20	21				
D2S1338	16	23				
D21S11	30	31.2				
D18S51	13	16				
D8S1179	11	17				

D3S1358	16	17				
D6S1043	18	19				
PENTAE	18	18				
D19S433	13	13				
PENTAD	9	13				
D1S1656	16	17.3				

4. Conclusion: The STR results of **TPC-1** cells are consistent with the genotypes of **TPC-1** cell lines in **EXPASY** database, the cell ID corresponded to **CVCL\_6298** , and the STR results **basically matched**. In the test, **No multi allele was found** in this cell line.

## 5. Attached Image

