

STR Authentication Report For Cell Line

1. Cell Name: U-87 MG

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: U-87 MG			Cell name: U-87MG		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	11	12		11	12	
D13S317	8	11		8	11	
D7S820	8	9		8	9	
D16S539	12	12		12	12	
VWA	15	17		15	17	
TH01	9.3	9.3		9.3	9.3	
AMEL	X	X		X	X	
TPOX	8	8		8	8	
CSF1PO	10	11		10	11	
D12S391	18	21				
FGA	18	24				
D2S1338	20	23				
D21S11	28	32.2				
D18S51	13	13				
D8S1179	10	11				

D3S1358	16	17				
D6S1043	11	18				
PENTAE	7	14				
D19S433	15	15.2				
PENTAD	9	14				
D1S1656	15	15				

4. Conclusion: The STR results of **U-87 MG** cells are consistent with the genotypes of **U-87 MG** cell lines in **DSMZ** database, the cell ID corresponded to **HTB-14** , and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

5. Attached Image



